

DOSI-FUSER®

Portable Elastomeric Pump

The DOSI-FUSER® is a single-use, continuous-infusion system for ambulatory patients that works without batteries or electricity. It consists of an elastomeric balloon inside a rigid, transparent container and an infusion line with a Luer Lock connector. As a result of the pressure from the elastomeric balloon, the medication circulates through the capillary element or the capillary tube which determines the flow rate regardless of the atmospheric pressure.



- Provides parenteral drug infusions at a constant flow without restricting patient mobility
- Suitable for ambulatory use
- May be used for intravenous, intra-arterial, epidural or subcutaneous administration
- Accuracy of infusion time is $\pm 15\%$
- Fluid level indicator
- 1.2 micron in-line particulate filter
- 0.02 micron air eliminating filter

Process Control. Protection. Peace of Mind.

Manufactured by: Leventon, S.A.U.