




High Frequency Sealer

- Seals PVC tubing on bags after filling providing a closed system
- No disposables needed to operate
- Standard power source

Printer

- Issues labels that identify content of prepared solution
- Documents results of filter integrity test
- Can print user-specified bar codes

Order No.	Description	Qty.
HFXS-01	High Frequency Sealer	—
PRT-01	Printer	—
LBL-12	Printer labels (small)	1000
LBL-14	Printer labels (large)	1000



Hanger Bags

- Provide simple, efficient way to hang vials of all sizes

Gri-fill® Cart

- Maximizes efficient utilization of space
- Moves system to unused space when not needed
- Organizes workflow and supplies

Order No.	Description	Qty.
HGR-06	Hanger stand (6-hook)	—
HBG-01	Hanger bags	1000
CRT-01	Gri-fill® Cart (37"w x 24"d x 32"h)	—