


The **Gri-fill®** System bags have the capability to filter every dose, offering increased confidence in delivery of a sterile product to your patients. The integral 0.22 micron filter on Gri-bag™, Gri-bag™ AP, and Gri-flex™ containers assures the maintenance of sterility in a closed system backed by an integrity test. The Gri-bag™ AP, with **Additive Port**, offers the flexibility to use multiple set spike sizes or to introduce another additive just prior to administration.

## Gri-bag™ and Gri-bag™ AP

- PVC bag with 0.22 micron filter in the filling port
- Maintains sterility of solution as it is filled through the port
- Cellulose acetate filter certified for retention of *Brevundimona (Pseudomonas) diminuta* at a minimum concentration of  $10^7$  organisms per  $cm^2$



**Gri-Bag AP** 500 ml  
Bags for Gri-Fill System

RS-LO—  
Sol A: MORPHINE Quan: 500mg Vol: 25.0ml  
Sol B: NaCl Quan: - Vol: 225.0ml  
Final container: Gri-bag Total volume: 250.0ml

LOT: CS1-80  
EXP12-05

User: SUPERV Date: 2005-11-22 14:32  
Integrity test: PASSED

LOT  
Laboratorios Grifols, S.A.  
Can Guasch, 2 - Parets del Vallès  
E08150 Barcelona - SPAIN 6004192  
**GRIFOLS**

Order No.	Description	Pack Factor	Case Qty.
GB100	Gri-bag™ 100mL	5 ea/unit	20
GB500	Gri-bag™ 500mL	5 ea/unit	20
GB1000	Gri-bag™ 1000mL	5 ea/unit	15
GBAP100	Gri-bag™ AP 100mL	5 ea/unit	20
GBAP500	Gri-bag™ AP 500mL	5 ea/unit	20
GBAP1000	Gri-bag™ AP 1000mL	5 ea/unit	15



## Gri-Flex™

- Polypropylene bag with 0.22 micron filter in the filling port
- Maintains sterility of solutions as it is filled through the port
- Supor™ polioethersulfona low protein-binding filter certified for retention of *Brevundimona (Pseudomonas) diminuta* at a minimum concentration of  $10^7$  organisms per  $\text{cm}^2$

Order No.	Description	Pack Factor	Case Qty.
GF100	Gri-flex™ 100mL	5 ea/unit	20
GF250	Gri-flex™ 250mL	5 ea/unit	20
GF500	Gri-flex™ 500mL	5 ea/unit	20
GF1000	Gri-flex™ 1000mL	5 ea/unit	15



Distributed in USA by:

**Grifols USA, Inc**

8880 NorthWest 18th Terrace  
Miami, FL 33172  
USA  
Tel.: [800] 379-0957  
[www.gri-fill.grifolsusa.com](http://www.gri-fill.grifolsusa.com)