

Optimize your use of the Gri-fill® System with MixManager.

The MixManager organizes your workflow and simplifies the order entry process while allowing you to link up to eight compounders for maximum efficiency. Reports help meet documentation requirements and streamline the workflow in compounding areas.

- Windows®-based
- Easy to use
- Provides data repository for compounded doses
- Helps organize compounding
- Links up to 8 **Gri-fill®** 3.0 compounders
- Creates templates to generate both batch and patient-specific doses
- Allows for customization of label output, solution characteristics and reconstitution mode

MixManager 1.0 Beta2

mixmanager GRIFOLS System

Order Entry - View Mode

Date: 08/12/2004 10:08:45

Patient ID:

Patient Name:

Bed Number: Location: TREATMENT 1 Cancer Care Center

Physician:

Doses Ordered: 1 Completed: 0 Cancelled: 0 Status: Verified

Drug Template	Drug	Dose	Volume
CHOP	Cyclophosphamide	1000 mg	10.5 ml
Prep.: T-21	0.9% Sodium Chloride		47.5 ml
No. of Doses: 1	Container: GR-BAG		Total vol: 100 ml

EXIT VIEW D. T. VIEW SCREEN VERIFY GO TO VER. STATUS PRINT BACK HELP

Order No.	Description
MMX4.0	MixManager Software

Distributed in USA by:

Grifols USA, Inc

8880 NorthWest 18th Terrace
Miami, FL 33172
USA
Tel.: [800] 379-0957
www.gri-fill.grifolsusa.com